


**EN** Titanium Forceps

**DE** Titan-Pinzetten

**ES** Pinzas de titanio

**Page 115 to 116**





**ATRAUMA**



**TITANIUM**



OI 10-2401-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 30 mm  
Jaw Length 9 mm  
Weight 1 gr



OI 10-2402-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 35 mm  
Jaw Length 13 mm  
Weight 1 gr



OI 10-2403-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 40 mm  
Jaw Length 14 mm  
Weight 1 gr



OI 10-2404-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 50 mm  
Jaw Length 23 mm  
Weight 1.5 gr



OI 10-2405-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 30 mm  
Jaw Length 9 mm  
Weight 1 gr



OI 10-2406-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 35 mm  
Jaw Length 13 mm  
Weight 1 gr



OI 10-2407-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria


Closing Pressure 220 +/- 10 gr  
Length 40 mm  
Jaw Length 14 mm  
Weight 1 gr



OI 10-2408-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 50 mm  
Jaw Length 23 mm  
Weight 1.5 gr



OI 10-2409-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 30 mm  
Jaw Length 9 mm  
Weight 1 gr



OI 10-2410-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 35 mm  
Jaw Length 13 mm  
Weight 1 gr



OI 10-2411-00  
Vena  
Closing pressure  
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr  
Length 40 mm  
Jaw Length 14 mm  
Weight 1 gr



OI 10-2412-00  
Vena  
Closing pressure  
110 +/- 10 gr


Arteria

Closing Pressure 220 +/- 10 gr  
Length 50 mm  
Jaw Length 23 mm  
Weight 1.5 gr






**ATRAUMA**



**DEBAKEY**



OI 10-2413-16  
16 cm/6 1/4"  
1.5 mm



OI 10-2413-20  
20.5 cm/8 1/8"  
1.5 mm



OI 10-2413-24  
24.5 cm/9 5/8"  
1.5 mm



OI 10-2414-16  
16 cm/6 1/4"  
1.8 mm



OI 10-2414-20  
20.5 cm/8 1/8"  
1.8 mm



OI 10-2414-24  
24.5 cm/9 5/8"  
1.8 mm

